# Your Liner<sup>®</sup>

Regulatory notes for professionals

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ID: 4.3.1.2.3\_Regulatory\_Notes\_CPOs\_SMT-L

# ENGLISH

1. Troubleshooting	3
2. Contraindications	3
3. Suspension type & compatible components	4
4. Liner reception & fitting	4
5. Miscellaneous	4
6. Ordering instructions	5

## 1. Troubleshooting

Problem	Likely cause	Possible solution
Seam is wearing or unraveling	Movement or rough socket edge	Apply instant adhesive to the seam
Silicone and fabric separate at the proximal edge	Liner is trimmed too close to the socket, and the socket edge rubs against the liner	Liner needs to be changed and trimmed longer
Fabric shows wear signs along socket trim line	The liner rubs against the trim line of the socket	Smooth out the edge of the trim line. Check that liner does not rub across the socket during gait. Optionally, add cover, like leather, to the socket edge.
Fabric or silicone is breaking down in the anterior distal tibia area	Socket has become loose in that area	Insert a pad in the socket in that area to tighten the fit.
Your patient experiences excessive sweating	Skin is not yet used to the closed environment that a liner imposes	The excessive sweating phenomenon should continue after a few weeks of using Your™ Liner. Patients that engage in physical activity that leads to sweating should periodically doff the liner and dry it off to prevent slippage of the liner on the limb. Some patients have tried to apply anti-perspirant to the residual limb before donning the liner, with effective results. Instructions shall be followed and anti-perspirant shall NOT be used if there are open sores on the skin. Please note that non alum-based products might damage the liner.
Patient experiences sweating, even though he has worn a Your™ Liner for a long time without excessive sweating	Some air manages to get into the liner	This may be due to limb shrinkage. Check the socket and liner fit. If the liner fit is loose, you need to order a new Your™ Liner with updated scan for new limb shape and volume. If the socket fit is loose, tighten up the socket.
	Liner presents a hole	Your™ Liner needs to be renewed
Top edge of the liner is rolling down the leg	Liner is too tight at the proximal end	The residual limb of your patient may have increased in volume. A new Your™ Liner shall be ordered with updated scan for new limb shape and volume.
	Liner is too short	The liner does not fit your patient's limb. There may have been an error during the scanning process. Contact MotionTech® customer service to check if the model provided matches your patient's limb measurements.

Problem	Likely cause	Possible solution
		The socket fit shall be tightened. If this is not efficient, a new Your™ Liner shall be ordered with updated scan for the new volume and shape of your patient's limb. To ensure that the irritation does not come from allergic reaction, apply a piece of silicone from a Your™ Liner in another area of the body and check for skin reaction.
Patient presents irritation all over the residual limb	Residual limb has shrunk and the liner or socket fit became loose	Ask your patient if he started using a new cleaning product, lotion, cream, powder or other skin conditioner. Hydrocarbon oils, animal fats and oils can often produce allergic reactions. Antibacterial soap has also been linked to skin irritations with some amputees.
	Use of an irritating soap lotion, cream, or other skin conditioner	
	Hygiene routine not properly implemented	Review the hygiene instructions with your patient, including cleaning, rinsing and disinfection procedures.
Patient presents irritation along the proximal edge of the liner	Patient does not follow the proper rolling procedure and pulls the liner up on the the skin instead.	Check your patient's donning technique. Reduced irritations have also been obtained by cutting the proximal edge of the liner in a wave pattern.
Open wounds or other skin injuries	Many possible causes. Socket fit and other specific patients conditions are most likely.	The specific cause need to be singled out. During healing, Your™ Liner should not be worn if the condition is too severe. If the injury is minor, bandages should be applied to cover the wound and prevent direct contact with Your™ Liner.
Residual limb and liner rotate inside the socket	Residual limb is highly cylindrical and socket is too loose	Tighten socket fit

## 2. Contraindications

### Contraindications for all Your® Liners

Although silicone liners have been used successfully with thousands of patients, MotionTech does not recommend using Your® Liners with the following types of patients:

- any patient who is either not capable or not willing to comply with the cleaning and hygiene procedures described in these Instructions for Use
- any patient who lacks the hand function to roll the liner successfully up onto the residual limb, and who does not have a caregiver to assist with this process
- any patient who lacks the ability to orient the liner properly on the limb, for example due to limited eyesight
- any patient who has a skin condition which could be aggravated by covering the limb with a non-porous liner, or who is known for developing allergic reactions to silicone

### Contraindications for Your® Tibial Suction Liners

Your® Tibial Suction Liners are intended to work safely when used with a dedicated prosthetic socket equipped with a standard expulsion valve. A Your® Tibial Suction Liner must also be used with a supplementary suspension system independent from the liner, in case the patient presents a short residual limb, substantial daily volume fluctuations, or if sufficient suspension can not be achieved without such supplementary suspension system. A Your® Tibial Suction Liner can not be provided without supplementary suspension method to a patient who presents a short residual limb or substantial residual limb volume fluctuations during the day, due to the risks of loss of suspension.

### Contraindications for Your® Femoral Suction Liners

Your® Femoral Suction Liners are intended to work safely when used with a dedicated prosthetic socket equipped with a standard expulsion valve. A Your® Femoral Suction Liner must also be used with a supplementary suspension system independent from the liner, in case the patient presents a short residual limb, substantial daily volume fluctuations, or if sufficient suspension can not be achieved without such supplementary suspension

system. A Your® Femoral Suction Liner can not be provided without supplementary suspension method to a patient who presents a short residual limb or substantial residual limb volume fluctuations during the day, due to the risks of loss of suspension.

#### **Contraindications for Your® Femoral Pin-lock Liners**

Your® Femoral Pin-lock Liners shall not be provided to patients with a very short residual limb (approximately 7cm (3") or shorter as measured from the perineum to the distal end along the medial side), unless a supplementary suspension system independent from the liner is used (for example a Silesian belt or a pelvic band), due to the risks of loss of suspension.

### **3. Suspension type & compatible components**

#### **Suspension type & components compatible with Your® Tibial Cushion Liners:**

A Your® Tibial Cushion Liner is designed to be used with a suspension sleeve. Even if you decide to add some anatomical suspension in your socket design - which is not recommended, as Your® Tibial Cushion Liners work best with total weight bearing sockets -, you MUST provide a suspension sleeve to your patient in order for the suspension to be efficient. Not using a suspension sleeve with a Your® Tibial Cushion Liners may result in slippage of the prosthesis on the liner.

It is your responsibility, as ordering prosthetist, to ensure that the suspension function provided by the chosen suspension sleeve is adequate and sufficient.

A Your® Tibial Cushion Liner does not hold any weight bearing or suspension function, and its strength is therefore not safety-related. Your® Tibial Cushion Liners have therefore not been subjected to strength tests. For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

#### **Suspension type & components compatible with Your® Tibial Suction Liners:**

A Your® Tibial Suction Liner is designed to be used with a dedicated socket equipped with a standard expulsion valve. Before choosing this type of liner for a patient, observe the contraindications presented in the present manual. It is mandatory that you ensure that the suspension strength obtained through the vacuum suspension is sufficient prior to providing the Your® Tibial Suction Liner to your patient. If sufficient suspension can not be reached through socket adaptations, you must provide him with a supplementary suspension system independent from the liner.

Extensive efforts have been implemented to ensure that the suspension strength of Your® Tibial Suction Liners is as high as possible. However, as Your® Tibial Suction Liners are tailor-made, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. Furthermore, the strength of the suspension depends to a great extent on the fit with the prosthetic socket. This is why the suspension function of Your® Tibial Suction Liners must be validated by the certified prosthetist prior to provision of such a device to the patient.

Your® Tibial Suction Liners have therefore only been subjected to shear fatigue tests to assess the strength of the suction lines. If any sign of excessive movement or product wear is detected on or around the suction lines of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

#### **Suspension type & components compatible with Your® Tibial Pin-lock Liners**

A Your® Tibial Pin-lock Liner is designed to be used with a distal attachment pin. MotionTech provides two options regarding compatible pin-lock systems: a M6 and a M10 thread. It is your responsibility, as ordering prosthetist, to ensure that the selected pin-lock system to be assembled on the prosthesis is compatible with the chosen connector (M6 or M10).

Apply Loctite 410/411 or equivalent compound to the threads of the pin and torque to 4 Nm (3 ft.lbs). It is not recommended to use a pin without a flange with Your® Tibial Pin-lock Liners. If you do, follow the same process but torque only to 1.5 Nm (1.1 ft.lbs).

Extensive efforts have been implemented to ensure that the suspension strength of Your® Tibial Pin-lock liners is as high as possible. Your® Tibial Pin-lock Liners have been tested under axial tensile load to assess their proof and ultimate strength under such loading conditions. However, as the shape of the locking system of Your® Tibial Pin-lock Liners are tailored to each patient, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. These conducted tests provided indications that a Your® Tibial Pin-lock Liner may be used with a prosthesis which weighs up to 10kg.

If any sign of excessive movement or product wear is detected on or around the locking part of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

#### **Suspension type & components compatible with Your® Femoral Suction Liners**

A Your® Femoral Suction Liner is designed to be used with a dedicated socket equipped with a standard expulsion valve. Before choosing this type of liner for a patient, observe the contraindications presented in the present manual. It is mandatory that you ensure that the suspension strength obtained through

the vacuum suspension is sufficient prior to providing the Your® Femoral Suction Liner to your patient. If sufficient suspension can not be reached through socket adaptations, you must provide him with a supplementary suspension system independent from the liner.

Extensive efforts have been implemented to ensure that the suspension strength of Your® Femoral Suction Liners is as high as possible. However, as Your® Femoral Suction Liners are tailor-made, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. Furthermore, the strength of the suspension depends to a great extent on the fit with the prosthetic socket. This is why the suspension function of Your® Femoral Suction Liners must be validated by the certified prosthetist prior to provision of such a device to the patient.

Your® Femoral Suction Liners have therefore only been subjected to shear fatigue tests to assess the strength of the suction lines. If any sign of excessive movement or product wear is detected on or around the suction lines of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

#### **Suspension type & components compatible with Your® Femoral Pin-lock Liners**

A Your® Femoral Pin-lock Liner is designed to be used with a distal attachment pin. MotionTech provides two options regarding compatible pin-lock systems: a M6 and a M10 thread. It is your responsibility, as ordering prosthetist, to ensure that the selected pin-lock system to be assembled on the prosthesis is compatible with the chosen connector (M6 or M10).

Apply Loctite 410/411 or equivalent compound to the threads of the pin and torque to 4 Nm (3 ft.lbs). It is not recommended to use a pin without a flange with Your® Femoral Pin-lock Liners. If you do, follow the same process but torque only to 1.5 Nm (1.1 ft.lbs).

Extensive efforts have been implemented to ensure that the suspension strength of Your® Femoral Pin-lock liners is as high as possible. Your® Femoral Pin-lock Liners have been tested under axial tensile load to assess their proof and ultimate strength under such loading conditions. However, as the shape of the locking system of Your® Femoral Pin-lock Liners are tailored to each patient, classical strength tests are not sufficient to ensure a flawless mechanical strength for each designed product. These conducted tests provided indications that a Your® Femoral Pin-lock Liner may be used with a prosthesis which weighs up to 10kg.

If any sign of excessive movement or product wear is detected on or around the locking part of the liner at any time, your patient should be instructed to immediately discontinue use of the device and consult you.

For further indications on maximum prosthesis weight and/or user weight, refer to the indications provided by the manufacturer of the components involved in the weight bearing and suspension functions.

### **4. Liner reception & fitting**

#### **Liner size check**

Despite all the quality checks we implement throughout our production, Your® Liners remain custom-made devices. It is therefore mandatory that you be present the first time Your® Liner is donned onto your patient's limb to ensure its suitability before providing it to your patient.

Pay attention to the observed tension of the liner once donned and make sure your patient does not feel unusual compression on his limb. Leave Your® Liner on your patient's limb for at least 5 minutes to ensure that the compression is neither too high nor too low. If there is a doubt on the compressive behaviour of Your® Liner, contact MotionTech's customer service prior to letting your patient go with Your® Liner. If air pockets or wrinkles persist to appear after five attempts of carefully donning and aligning the liner, the liner should be deemed not adapted, and you should contact MotionTech' customer service for further indications. Similarly, if the liner is too tight, too loose or too short, you should contact MotionTech' customer service to assess the situation.

#### **Trimming the liner**

You may trim Your® Liner to a shorter length. For a BK liner, do not trim it shorter than 8cm (3.1") above the proximal end of the socket. Use sharp scissors, a paper cutter, or a dedicated Hand-Held Cutting Wheel. To help prevent the seam from unraveling: Place a small drop of instant adhesive on the seam directly below the trimmed edge.

## **5. Miscellaneous**

#### **General Disclaimer**

A Your® Liner can be order by a certified prosthetist only.

A Your® Liner offers the possibility to the prosthetist to customise a wide variety of parameters, in order to provide optimal comfort and performances to his-her patient. It is the responsibility of the ordering prosthetist to select the suitable parameters depending on the patient he intends to provide a Your® Liner with.

The Your® Liner documentation gives advices on which parameters should be selected for particular situations. They should be treated as guidelines, and not absolute rules, as each patient may differ. MotionTech's customer support teams are entirely at your disposal for further guidance on the liner design process.

MotionTech retains all intellectual property rights reflected in its physical products, regardless of the transfer of the physical products to another party or parties.

**Liability**

Swiss Motion Technologies SA, as manufacturer, recommends using Your® Liners only under the specified conditions and for the intended purposes. Your® Liners must be maintained according to the instructions for use. The manufacturer is not liable for any damage caused by combinations of component that were not specifically recommended by the manufacturer in the present document.

**6. Ordering instructions**

Important instructions for placing correct orders is provided in SMT's Frequently Asked Questions for professionals page (<https://www.motiontech.ch/faq-pros>).

SMT's ordering platform also include instructions and guidance on important steps to follow when placing orders (<https://www.motiontech.ch/yourlinerorderingplatform>).

Failure to properly follow these instructions may result in exposing your patients to risk of skin injuries such as blisters or other types of skin breakdown. Contact your local SMT representative for further guidance on proper processes for ordering Your® Liner.

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